

FOREST

control heads



IC CONTROL HEAD

- Ultra reliable pneumatic master control block. A thick cap prevents warping. Independent water relay valves remove water from the block eliminating destructive deposits. The result is a trouble-free life.
- Ultra fine water and air coolant adjustments.
- Hinged lid for easy access during routine maintenance.
- 1/2" holder bar can be lengthened to accept numerous ancillary products.
- Ergonomically placed brake release handle for precise positioning of control head. Euro control head has dual brake release for left and right use.
- Hygienist head conveniently incorporates HVE, saliva ejector and swing-away vacuum canister with solids collector.
- Optional designer friendly accents are available to compliment your sense of style.

CONTOURED IC



STANDARD IC
+ HYGIENIST IC



EURO




World Class Design : Beautiful | Functional | Affordable

Details on reverse




forest
AFFORDABLE EXCELLENCE




CONTOURED IC

- Plastic molded hinged lid and chassis
- Ultra-reliable pneumatic master control block
- Ultra-fine individual water and air coolant adjustment
- Push button flex arm brake release
- Pivot has dual brake release for left/right use



STANDARD & HYGIENIST IC

- Powder coated hinged metal lid
- Plastic molded chassis
- Ultra-reliable pneumatic master control block
- Ultra-fine individual water and air coolant adjustment
- Push button flex arm brake release
- Pivot has dual brake release for left/right use



EURO

- Plastic molded cover and handpiece pad
- Arm mounted tray holder
- Ultra-reliable pneumatic master control block
- Ultra-fine individual water and air coolant adjustment
- Dual push button flex arm brake release for left/right use

Mounting Option	Contoured Model/Price		Standard Model/Price		Standard Hyg. Model/Price		Euro Model/Price	
	Chair	8350	\$ 5965 *	8300	\$ 5965 *	4455	\$ 5865 *	8800
Cart	8138	\$ 3970	4338	\$ 3970	4548	\$ 4470	4165	\$ 4925
	4585	\$ 4810	4485	\$ 4810	4495	\$ 5210	4985	\$ 5765
Wall/Cabinet	3535	\$ 3720	3235	\$ 3720	3245	\$ 4220		
	3530	\$ 3720	3230	\$ 3720	3240	\$ 4220		
near under	2450	\$ 2380	2400	\$ 2380				
	3650	\$ 5775	3600	\$ 5775				
	5850	\$ 4955	5800	\$ 4955				

* includes sidebox

1.800.423.3555
www.forestmed.com

